

SUPPLY ORDER FORM

Facility Name:	Facility ID:
Facility Address:	Attention:
City, State, Zip:	Email:
Phone Number:	Order Date:
_	

FAX ORDER FORMS BACK TO (517)787-3801

SUPPLIES ARE GENERALLY SHIPPED VIA UPS. PLEASE ALLOW 5 BUSINESS DAYS FOR DELIVERY INDICATE THE QUANTITY OF EACH ITEM IN THE APPROPRIATE COLUMN BELOW

BLOOD COLLECTION TUBES						
Ind Qty	Pkg Qty	Product Description	Qty/Pkg			
		Red/Grey Top SST 8.5ml	100			
		Lavender Top Tube 4ml (K ² EDTA)	100			
		Light Blue Top Sodium Citrate 4.5ml	100			
		Red Top Plain 6ml				
		PPT (Plasma Preperation Tube) 8ml				
	_	Gray Top 4ml	_			
	_	Royal Blue Top EDTA 7ml				
	 Royal Blue Top - No Additive 6ml 					
	Yellow Top ACD "A" Solution 8.5ml					
	_	Green Top Sodium Heparin 10ml	_			
_		Blood Culture Set	2			
URINE COLLECTION CONTAINERS						
	-	24 Hr. Rigid - Plastic	_			
	_	60ml Sterile Cup	_			
	_	5ml Urine Transport - Gray Top	_			
OTHER SPECIMEN COLLECTION						
		Specimen Bag - Small	100			
		Specimen Bag - Large	100			
		Specimen Transport Boxes w/ Foam	24			
	_	Transport Tubes - Clear Plastic	_			
		UPS Shipping Bags - Laboratory Pak	50			
	_	UPS Mailing Labels	_			
NEEDLES						
_		21G X 1 1/4	48			
-		22G X 1 1/4	48			
-		23G X 3/4 Butterfly	50			
_		21G Butterfly	50			

	TRANSPORT MEDIA
Ind Qty	Product Description
	M4 Multi-Microbe Transport Kit
	Ova & Parasite Transport Vials
	C&S Transport - Stool
	Trichosel Broth
	Formalin Cont Small Large (circle one)
	OTHER SPECIMEN COLLECTION
	Culturette - Anaerobic
	Culturette - Aerobic (Yellow Top/Wound)
	Gen-Probe® Collection - Unisex
	Thin Prep - Pap Test
	Traditional Pap Kits - Slides
	Glucola - 50 75 100 (circle one)
	Stool for Occult Blood Card
	Culturette - Flu
	FORMS AND REQUISITIONS
	Client Supply Order Form
	Requisition Forms - Standard
	Cyto/Pathology Forms
	Blank Requisition Paper w/ Sticker Labels
	OTHER SUPPLIES
	Single Use Needle Holder - Portex
	Amber Tube 4ml - Plastic
	Tourniquets - Standard
·	Tourniquets - Latex Free
	Alcohol Prep Pads
	Autoclave Spore Kit
·	Cotton Balls

Additional Comments:		